


## Product Datasheet

### Human SIK2 ELISA Kit (orb563657)

<b>Description</b>	Human SIK2 ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	0.312-20ng/ml
<b>Target</b>	SIK2
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Application notes</b>	Competitive ELISA, Coated with Antigen
<b>Uniprot ID</b>	<b>Q9H0K1</b>
<b>Sensitivity</b>	0.188ng/ml
<b>Expiration Date</b>	Please enquire.



Human  
SIK2(Serine/threonine-  
protein kina...