


Product Datasheet

Human RSAD2 ELISA Kit (orb563671)

Description	Human RSAD2 ELISA Kit
Reactivity	Human
Range	0.78-50ng/ml
Target	RSAD2
Storage	2-8 °C for 6 months
Note	For research use only
Application notes	Competitive ELISA, Coated with Antigen
Uniprot ID	Q8WXG1
Sensitivity	0.469ng/ml
Expiration Date	Please enquire.



Human
RSAD2(Radical
S-adenosyl
methionin...