


## Product Datasheet

### Human CYP2C19 ELISA Kit (orb563717)

<b>Description</b>	Human CYP2C19 ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	0.156-10ng/ml
<b>Target</b>	CYP2C19
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Sample Types</b>	Serum,Plasma,Cell Culture Supernatant,Cell or tissue lysate,Other liquid samples
<b>Uniprot ID</b>	<b>P33261</b>
<b>Sensitivity</b>	0.094ng/ml
<b>Expiration Date</b>	Please enquire.



Human  
CYP2C19(Cytochrome  
P450 2C19) ELIS...