


## Product Datasheet

### Human BIRC7 ELISA Kit (orb563978)

|                        |                        |
|------------------------|------------------------|
| <b>Description</b>     | Human BIRC7 ELISA Kit  |
| <b>Reactivity</b>      | Human                  |
| <b>Range</b>           | 125-8000pg/ml          |
| <b>Target</b>          | BIRC7                  |
| <b>Storage</b>         | 2-8 °C for 6 months    |
| <b>Note</b>            | For research use only  |
| <b>Uniprot ID</b>      | <a href="#">Q96CA5</a> |
| <b>Sensitivity</b>     | 75pg/ml                |
| <b>Expiration Date</b> | Please enquire.        |



Human  
BIRC7(Baculoviral  
IAP repeat-  
conta...