
Product Datasheet

beta Arrestin 1(phospho-Ser412) antibody (orb5640)

Description

beta Arrestin 1(phospho-Ser412)

Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human beta-Arrestin 1 around the phosphorylation site of Ser412 (TG(p-S)PR)
Target	Phospho-beta-Arrestin 1 (Ser412)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	45 kDa
Uniprot ID	P49407
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.