


## Product Datasheet

### Human COL15 ELISA Kit (orb564059)

<b>Description</b>	Human COL15 ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	0.313-20ng/ml
<b>Target</b>	COL15
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Sensitivity</b>	0.188ng/ml
<b>Expiration Date</b>	Please enquire.



Human  
COL15(Collagen  
Type XV) ELISA  
Kit