


## Product Datasheet

### Human ERK ELISA Kit (orb564132)

<b>Description</b>	Human ERK(ELK-Related Tyrosine Kinase) ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	0.625-40ng/ml
<b>Target</b>	ERK
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>P29323</b>
<b>Sensitivity</b>	0.375ng/ml
<b>Expiration Date</b>	Please enquire.



Human ERK  
(ELK-  
Related  
Tyrosine  
Kinase) ...