


Product Datasheet

Human Flt3L ELISA Kit (orb564176)

Description	Human Flt3L ELISA Kit
Reactivity	Human
Range	15.625-1000pg/ml
Target	Flt3L
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum,Plasma,Cell Culture Supernatant,Cell or tissue lysate,Other liquid samples
Uniprot ID	P49771
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Human
Flt3L(FMS
Like Tyrosine
Kinase 3 L...