


Product Datasheet

Human ITG α V β 3(Integrin α V β 3) ELISA Kit (orb564178)

Description	Human ITG α V β 3(Integrin α V β 3) ELISA Kit
Reactivity	Human
Range	31.2-2000pg/ml
Target	ITG α V β 3
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	P06756, P05106
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Human ITG
 α V β 3(Integrin
 α V β 3) EL...