


Product Datasheet

Human fT3 ELISA Kit (orb564184)

| | |
|--------------------------|--|
| Description | Human fT3 ELISA Kit |
| Reactivity | Human |
| Range | 0.781-50pg/ml |
| Target | fT3 |
| Storage | 2-8 °C for 6 months |
| Note | For research use only |
| Application notes | Competitive ELISA, Coated with Antigen |
| Sensitivity | 0.469pg/ml |
| Expiration Date | Please enquire. |



Human fT3(Free
Triiodothyronine)
ELISA K...