


## Product Datasheet

### Human F-TESTO ELISA Kit (orb564186)

<b>Description</b>	Human F-TESTO ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	15.625-1000pg/ml
<b>Target</b>	F-TESTO
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Application notes</b>	Competitive ELISA, Coated with Antigen
<b>Sensitivity</b>	9.375pg/ml
<b>Expiration Date</b>	Please enquire.



Human F-TESTO(Free Testosterone) ELISA Ki...