


## Product Datasheet

### Human GAL8 ELISA Kit (orb564193)

<b>Description</b>	Human GAL8 ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	0.313-20ng/ml
<b>Target</b>	GAL8
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>000214</b>
<b>Sensitivity</b>	0.188ng/ml
<b>Expiration Date</b>	Please enquire.



Human  
GAL8(Galectin  
8) ELISA Kit