


## Product Datasheet

### Human GnRH ELISA Kit (orb564216)

<b>Description</b>	Human GnRH ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	15.625-1000pg/ml
<b>Target</b>	GnRH
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Application notes</b>	Competitive ELISA, Coated with Antigen
<b>Uniprot ID</b>	<b>P01148</b>
<b>Sensitivity</b>	9.375pg/ml
<b>Expiration Date</b>	Please enquire.



Human  
GnRH(Gonadotropin  
Releasing  
Hormon...