


Product Datasheet

Human Hepc ELISA Kit (orb564254)

| | |
|------------------------|-----------------------|
| Description | Human Hepc ELISA Kit |
| Reactivity | Human |
| Range | 31.25-2000pg/ml |
| Target | Hepc |
| Storage | 2-8 °C for 6 months |
| Note | For research use only |
| Uniprot ID | P81172 |
| Sensitivity | 18.75 pg/ml |
| Expiration Date | Please enquire. |



Human
Hepc(Hepcidin)
ELISA Kit