


## Product Datasheet

### Human IGF2BP3 ELISA Kit (orb564275)

<b>Description</b>	Human IGF2BP3 ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	0.313-20ng/ml
<b>Target</b>	IGF2BP3
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>000425</b>
<b>Sensitivity</b>	0.188ng/ml
<b>Expiration Date</b>	Please enquire.



Human  
IGF2BP3(Insulin  
Like Growth  
Factor...