


## Product Datasheet

### Human IGF2R ELISA Kit (orb564276)

<b>Description</b>	Human IGF2R ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	0.625-40ng/ml
<b>Target</b>	IGF2R
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<a href="#">P11717</a>
<b>Sensitivity</b>	0.375ng/ml
<b>Expiration Date</b>	Please enquire.



Human  
IGF2R(Insulin  
Like Growth  
Factor 2...