


## Product Datasheet

### Human IL-17 ELISA Kit (orb564284)

|                        |  |
|------------------------|--|
| <b>Description</b>     | Human IL-17A(Interleukin 17 A) ELISA Kit |
| <b>Reactivity</b>      | Human                                    |
| <b>Range</b>           | 31.25-2000pg/ml                          |
| <b>Target</b>          | IL-17                                    |
| <b>Storage</b>         | 2-8 °C for 6 months                      |
| <b>Note</b>            | For research use only                    |
| <b>Uniprot ID</b>      | <b>Q16552</b>                            |
| <b>Sensitivity</b>     | 18.75pg/ml                               |
| <b>Expiration Date</b> | Please enquire.                          |



Human IL-17(Interleukin 17) ELISA Kit