


## Product Datasheet

### Human LDHB ELISA Kit (orb564298)

<b>Description</b>	Human LDHB ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	3.125-200ng/ml
<b>Target</b>	LDHB
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>P07195</b>
<b>Sensitivity</b>	1.875ng/ml
<b>Expiration Date</b>	Please enquire.



Human  
LDHB(Lactate  
Dehydrogenase  
B) ELIS...