


## Product Datasheet

### Human MIP4alpha ELISA Kit (orb564324)

<b>Description</b>	Human MIP4alpha(Macrophage Inflammatory Protein 4 Alpha) ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	31.25-2000pg/ml
<b>Target</b>	MIP4 $\alpha$
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>Q9Y258</b>
<b>Sensitivity</b>	18.75pg/ml
<b>Expiration Date</b>	Please enquire.



Human  
MIP4 $\alpha$ (Macrophage  
Inflammatory  
Pro...