


## Product Datasheet

### Human N-MID-OT ELISA Kit (orb564356)

<b>Description</b>	Human N-MID-OT ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	0.625-40ng/ml
<b>Target</b>	N-MID-OT
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>P02818</b>
<b>Sensitivity</b>	0.375ng/ml
<b>Expiration Date</b>	Please enquire.



Human N-MID-OT(N-MID Osteocalcin) ELISA ...