


## Product Datasheet

# Human PIIIIP(Procollagen III C-Terminal ProPeptide) ELISA Kit (orb564380)

<b>Description</b>	Human PIIIIP(Procollagen III C-Terminal ProPeptide) ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	12.5-800pg/ml
<b>Target</b>	PIIIP
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Sensitivity</b>	7.5pg/ml
<b>Expiration Date</b>	Please enquire.



Human PIIIIP(Procollagen III C-Terminal ...