


## Product Datasheet

### Human PGRP-S ELISA Kit (orb564422)

<b>Description</b>	Human PGRP-S ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	0.313-20ng/ml
<b>Target</b>	PGRP-S
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>075594</b>
<b>Sensitivity</b>	0.188ng/ml
<b>Expiration Date</b>	Please enquire.



Human PGRP-S (Peptidoglycan Recognition P...