


## Product Datasheet

### Human PPP1R1A ELISA Kit (orb564459)

<b>Description</b>	Human PPP1R1A ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	62.5-4000pg/ml
<b>Target</b>	PPP1R1A
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>Q13522</b>
<b>Sensitivity</b>	37.5pg/ml
<b>Expiration Date</b>	Please enquire.



Human  
PPP1R1A(Protein  
Phosphatase 1  
Regu...