


Product Datasheet

Human SPD ELISA Kit (orb564537)

Description	Human SPD ELISA Kit
Reactivity	Human
Range	1.563-100ng/ml
Target	SPD
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum,Plasma,Cell Culture Supernatant,Cell or tissue lysate,Other liquid samples
Uniprot ID	P35247
Sensitivity	0.938ng/ml
Expiration Date	Please enquire.



Human
SPD(Pulmonary
Surfactant
Associate...