


Product Datasheet

Human sTREM-1 ELISA Kit (orb564548)

| | |
|------------------------|-------------------------|
| Description | Human sTREM-1 ELISA Kit |
| Reactivity | Human |
| Range | 31.25-2000pg/ml |
| Target | sTREM-1 |
| Storage | 2-8 °C for 6 months |
| Note | For research use only |
| Uniprot ID | Q9NP99 |
| Sensitivity | 18.75pg/ml |
| Expiration Date | Please enquire. |



Human
sTREM-
1(soluble
Triggering
Recepto...