


## Product Datasheet

### Human VEGF-D ELISA Kit (orb564603)

<b>Description</b>	Human VEGF-D ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	62.5-4000pg/ml
<b>Target</b>	VEGF-D
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>O43915</b>
<b>Sensitivity</b>	37.5pg/ml
<b>Expiration Date</b>	Please enquire.



Human  
VEGF-  
D(Vascular  
Endothelial  
Growth...