


## Product Datasheet

### Human alpha2-AP ELISA Kit (orb564607)

<b>Description</b>	Human alpha2-AP(alpha2-Antiplasmin) ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	12.5-800ng/ml
<b>Target</b>	$\alpha$ 2-AP
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>P08697</b>
<b>Sensitivity</b>	7.5ng/ml
<b>Expiration Date</b>	Please enquire.



Human  $\alpha$ 2-  
AP( $\alpha$ 2-  
Antiplasmin)  
ELISA Kit