


Product Datasheet

Human OT ELISA Kit (orb564761)

Description	Human OT ELISA Kit
Reactivity	Human
Range	15.625-1000pg/ml
Target	OT
Storage	2-8 °C for 6 months
Note	For research use only
Application notes	Competitive ELISA, Coated with Antigen
Sample Types	Serum,Plasma,Cell Culture Supernatant,Cell or tissue lysate,Other liquid samples
Uniprot ID	P01178
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Human
OT(Oxytocin)
ELISA Kit