


Product Datasheet

Human SP ELISA Kit (orb564786)

Description	Human SP ELISA Kit
Reactivity	Human
Range	7.8-500pg/ml
Target	SP
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum,Plasma,Cell Culture Supernatant,Cell or tissue lysate,Other liquid samples
Uniprot ID	P20366
Sensitivity	4.69pg/ml
Expiration Date	Please enquire.



Human
SP(Substance
P) ELISA Kit