


## Product Datasheet

### Human ITGA1 ELISA Kit (orb564793)

<b>Description</b>	Human ITGA1 ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	0.313-20ng/ml
<b>Target</b>	ITGA1/CD49a
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>P56199</b>
<b>Sensitivity</b>	0.188ng/ml
<b>Expiration Date</b>	Please enquire.



Human  
ITGA1/CD49a(Integrin  
alpha-1) ELIS...