


Product Datasheet

Human NEXN ELISA Kit (orb564859)

Description	Human NEXN ELISA Kit
Reactivity	Human
Range	31.25-2000pg/ml
Target	NEXN
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum,Plasma,Cell Culture Supernatant,Cell or tissue lysate,Other liquid samples
Uniprot ID	Q0ZGT2
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Human
NEXN(Nexilin)
ELISA Kit