


Product Datasheet

Human ARNTL ELISA Kit (orb564942)

Description	Human ARNTL (Aryl hydrocarbon receptor nuclear translocator-like protein 1) ELISA Kit
Reactivity	Human
Range	15.625-1000pg/ml
Target	ARNTL
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	000327
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Human
ARNTL(Aryl
hydrocarbon
receptor
nu...