


## Product Datasheet

### Human HRP 2 ELISA Kit (orb565007)

<b>Description</b>	Human HRP 2 ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	31.2-2000pg/ml
<b>Target</b>	HRP 2
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<a href="#">Q7Z4V5</a>
<b>Sensitivity</b>	18.75pg/ml
<b>Expiration Date</b>	Please enquire.



Human HRP 2 (Hepatoma-derived growth factor...)