


Product Datasheet

Mouse NT-3 ELISA Kit (orb565076)

Description	Mouse NT-3 ELISA Kit
Reactivity	Mouse
Range	31.2-2000pg/ml
Target	NT-3
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum,Plasma,Cell Culture Supernatant,Cell or tissue lysate,Other liquid samples
Uniprot ID	P20181
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Mouse NT-3 (Neurotrophin 3) ELISA Kit