


## Product Datasheet

### Mouse Angiopoietin 2 ELISA Kit (orb565080)

|                        |                                |
|------------------------|--------------------------------|
| <b>Description</b>     | Mouse Angiopoietin 2 ELISA Kit |
| <b>Reactivity</b>      | Mouse                          |
| <b>Range</b>           | 39.1-2500pg/ml                 |
| <b>Target</b>          | Ang2                           |
| <b>Storage</b>         | 2-8 °C for 6 months            |
| <b>Note</b>            | For research use only          |
| <b>Uniprot ID</b>      | <b>Q64438</b>                  |
| <b>Sensitivity</b>     | 22.5pg/ml                      |
| <b>Expiration Date</b> | Please enquire.                |



Mouse  
Ang2(Angiopoietin-  
2) ELISA Kit