


## Product Datasheet

### Mouse Flt3L ELISA Kit (orb565106)

|                          |                                                                                                                                                       |
|--------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>       | Mouse Flt3L ELISA Kit                                                                                                                                 |
| <b>Reactivity</b>        | Mouse                                                                                                                                                 |
| <b>Range</b>             | 15.625-1000pg/ml                                                                                                                                      |
| <b>Target</b>            | Flt3L                                                                                                                                                 |
| <b>Storage</b>           | 2-8 °C for 6 months                                                                                                                                   |
| <b>Note</b>              | For research use only                                                                                                                                 |
| <b>Application notes</b> | Flt3L ELISA Kit allows for the in vitro quantitative determination of Flt3L concentrations in plasma, tissue homogenates and other biological fluids. |
| <b>Uniprot ID</b>        | <b>P49772</b>                                                                                                                                         |
| <b>Sensitivity</b>       | 9.375pg/ml                                                                                                                                            |
| <b>Expiration Date</b>   | Please enquire.                                                                                                                                       |



Mouse Flt3L  
ELISA Kit