


## Product Datasheet

### Mouse PRGN ELISA Kit (orb565151)

<b>Description</b>	Mouse PRGN ELISA Kit
<b>Reactivity</b>	Mouse
<b>Range</b>	0.078-5ng/ml
<b>Target</b>	PGRN
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>P28798</b>
<b>Sensitivity</b>	0.047ng/ml
<b>Expiration Date</b>	Please enquire.



Mouse  
PRGN(Progranulin)  
ELISA Kit