


## Product Datasheet

### Mouse RBP4 ELISA Kit (orb565157)

<b>Description</b>	Mouse RBP4 ELISA Kit
<b>Reactivity</b>	Mouse
<b>Range</b>	1.563-100ng/ml
<b>Target</b>	RBP4
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>Q00724</b>
<b>Sensitivity</b>	0.938ng/ml
<b>Expiration Date</b>	Please enquire.



Mouse  
RBP4(Retinol  
Binding  
Protein 4,  
Pl...