


Product Datasheet

Mouse TLR2 ELISA Kit (orb565169)

Description	Mouse TLR2 ELISA Kit
Reactivity	Mouse
Range	0.312-20ng/ml
Target	TLR2
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum,Plasma,Cell Culture Supernatant,Cell or tissue lysate,Other liquid samples
Uniprot ID	Q9QUN7
Sensitivity	0.188ng/ml
Expiration Date	Please enquire.



Mouse
TLR2(Toll-like
receptor 2)
ELISA K...