


## Product Datasheet

### Mouse Inhba ELISA Kit (orb565214)

<b>Description</b>	Mouse Inhba ELISA Kit
<b>Reactivity</b>	Mouse
<b>Range</b>	15.625-1000pg/ml
<b>Target</b>	Inhba
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>Q04998</b>
<b>Sensitivity</b>	9.375pg/ml
<b>Expiration Date</b>	Please enquire.



Mouse  
Inhba(Inhibin  
beta A  
chain) ELISA  
...