


Product Datasheet

Mouse Endothelin 1 ELISA Kit (orb565244)

Description	Mouse Endothelin 1 ELISA Kit
Reactivity	Mouse
Range	1.25-80pg/ml
Target	Edn1
Storage	2-8 °C for 6 months
Note	For research use only
Application notes	Competitive ELISA, Coated with Antigen
Uniprot ID	P22387
Sensitivity	0.625pg/ml
Expiration Date	Please enquire.



Mouse
Edn1(Endothelin-
1) ELISA Kit