


## Product Datasheet

### Mouse Gpx1 ELISA Kit (orb565279)

<b>Description</b>	Mouse Gpx1 ELISA Kit
<b>Reactivity</b>	Mouse
<b>Range</b>	4.69-300IU/ml
<b>Target</b>	Gpx1
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>P11352</b>
<b>Sensitivity</b>	2.813IU/ml
<b>Expiration Date</b>	Please enquire.



Mouse  
Gpx1(Glutathione  
peroxidase 1)  
ELI...