


## Product Datasheet

### Mouse Hmgb1 ELISA Kit (orb565293)

|                        |                       |
|------------------------|-----------------------|
| <b>Description</b>     | Mouse Hmgb1 ELISA Kit |
| <b>Reactivity</b>      | Mouse                 |
| <b>Range</b>           | 15.625-1000pg/ml      |
| <b>Target</b>          | Hmgb1                 |
| <b>Storage</b>         | 2-8 °C for 6 months   |
| <b>Note</b>            | For research use only |
| <b>Uniprot ID</b>      | <b>P63158</b>         |
| <b>Sensitivity</b>     | 9.375pg/ml            |
| <b>Expiration Date</b> | Please enquire.       |



Mouse  
Hmgb1(High  
mobility  
group  
protein ...