


## Product Datasheet

### Mouse ADIPOR1 ELISA Kit (orb565294)

<b>Description</b>	Mouse ADIPOR1 ELISA Kit
<b>Reactivity</b>	Mouse
<b>Range</b>	3.125-200ng/ml
<b>Target</b>	ADIPOR1
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>Q91VH1</b>
<b>Sensitivity</b>	1.875ng/ml
<b>Expiration Date</b>	Please enquire.



Mouse  
ADIPOR1(Adiponectin  
Receptor 1) EL...