


## Product Datasheet

### Mouse CYP27B1 ELISA Kit (orb565304)

<b>Description</b>	Mouse CYP27B1 ELISA Kit
<b>Reactivity</b>	Mouse
<b>Range</b>	0.312-20ng/ml
<b>Target</b>	Vdr
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Sample Types</b>	Cell or tissue lysate,Other liquid samples
<b>Uniprot ID</b>	<b>P48281</b>
<b>Sensitivity</b>	0.188ng/ml
<b>Expiration Date</b>	Please enquire.



Mouse  
Vdr(Vitamin  
D3  
receptor)  
ELISA Kit