


## Product Datasheet

### Mouse CD274 ELISA Kit (orb565308)

<b>Description</b>	Mouse CD274 ELISA Kit
<b>Reactivity</b>	Mouse
<b>Range</b>	31.2-2000pg/ml
<b>Target</b>	CD274
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>Q9EP73</b>
<b>Sensitivity</b>	18.75pg/ml
<b>Expiration Date</b>	Please enquire.



Mouse  
CD274(Programmed  
cell death 1 liga...