


## Product Datasheet

### Mouse ACTH ELISA Kit (orb565358)

<b>Description</b>	Mouse ACTH ELISA Kit
<b>Reactivity</b>	Mouse
<b>Range</b>	12.5-800pg/ml
<b>Target</b>	ACTH
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Application notes</b>	Competitive ELISA, Coated with Antigen
<b>Uniprot ID</b>	<b>P01193</b>
<b>Sensitivity</b>	7.5pg/ml
<b>Expiration Date</b>	Please enquire.



Mouse  
ACTH(Corticotropin)  
ELISA Kit