


## Product Datasheet

### Mouse Lyve1 ELISA Kit (orb565463)

<b>Description</b>	Mouse Lyve1 ELISA Kit
<b>Reactivity</b>	Mouse
<b>Range</b>	78-5000pg/ml
<b>Target</b>	Lyve1
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>Q8BHC0</b>
<b>Sensitivity</b>	46.9pg/ml
<b>Expiration Date</b>	Please enquire.



Mouse  
Lyve1(Lymphatic  
vessel  
endothelial...