


## Product Datasheet

### Mouse Acly ELISA Kit (orb565473)

<b>Description</b>	Mouse Acly ELISA Kit
<b>Reactivity</b>	Mouse
<b>Range</b>	0.312-20ng/ml
<b>Target</b>	Acly
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<a href="#">Q91V92</a>
<b>Sensitivity</b>	0.188ng/ml
<b>Expiration Date</b>	Please enquire.



Mouse  
Acly(ATP-  
citrate  
synthase)  
ELISA K...