


Product Datasheet

Mouse Adam28 ELISA Kit (orb565483)

| | |
|------------------------|------------------------|
| Description | Mouse Adam28 ELISA Kit |
| Reactivity | Mouse |
| Range | 31.2-2000pg/ml |
| Target | Adam28 |
| Storage | 2-8 °C for 6 months |
| Note | For research use only |
| Uniprot ID | Q9JLN6 |
| Sensitivity | 18.75pg/ml |
| Expiration Date | Please enquire. |



Mouse
Adam28(Disintegrin
and metalloprot...